

from the European Synchrotron Radiation Facility (ESRF) in Grenoble and Hyotcherl Ihee from Korea's Institute for Basic Science and the Korea Advanced Institute of Science and Technology (KAIST).

The first moments of a chemical reaction


Molecules are usually born in high energy states which are accompanied by large amplitude vibrations in the bonds between the atoms. In solution, a newly formed molecule cools down by transferring energy to the surrounding solvent molecules. This process takes place in just a few picoseconds - a few trillionths of a second - and has so far

never been directly visualised. Solvent plays an important role in solution-phase chemical reactions by serving as an energy source to activate the reaction as well as a heat bath to stabilise the products. As a result, the properties of the solvent significantly affect the energy landscape, rates and pathways of a reaction in solution. How solute and solvent molecules interact and the effect of this interaction on the outcome of chemical reactions has been a topic of intense research in the field of reaction dynamics for several decades. It is still challenging experimentally to probe the bonds that are formed during a chemical reaction, especially the changes between solute and solvent accompanying the structural change of reacting solute molecules. The challenge for these experiments is to develop new tools that have both sufficient time resolution and structural sensitivity.

Unprecedented time resolution

The novel technique of time slicing developed by the researchers consists in using a very short 1 picosecond laser pulse, a hundred times shorter than the X-ray pulse. The laser pulse dissociates the molecule into two atoms that are trapped in the solvent cage. Then a delayed X-ray pulse probes their reassociation at that time delay. The delay is controlled to a precision of 1 ps. The association process is probed in10 ps time steps and the data treated with a new deconvolution algorithm to sharpen the data.

By combining this technique with time-resolved X-ray scattering, the team was able to reach an unprecedented time resolution of 10 ps. Coupled with structural sensitivity offered by the very high pulse intensity from a quasi monochromatic undulator (U17) and a very efficient detector (FReLoN), the scientists were able to measure the dynamics of atom dissociation and reassociation and the vibrational relaxation of 12 in two different solvents.



Researchers filmed the dissociation of the iodine molecule and the subsequent reassociation and vibrational relaxation of the molecule.



02-04-2012

more






E Like K 0



The experimental setup showing schematically the laser and x-ray pulse in space 100 picosecond before they strike the sample. The 100 ps x-ray pulse forms a 30 mm long needle with a 0.3 mm long laser pulse inside it. By varying the relative position of the pulses and using deconvolution algorithms, the distance between the iodine atoms can be measured with 10 ps resolution.

The team found that the I-I bond distance increases from 2.65 Å to a maximum of about 4 Å before the iodine atoms draw back together. The distance differs depending on the solvent: it is greater in cyclohexane than in carbon tetrachloride, which has a higher molecular weight.

A film of the vibrating atoms inside the solvent cage was finally made by stitching together 30 snapshots taken at delays starting from 0 ps to 300 ps. Click on the photo below to view the animation.



The animation shows the dissociation and association of iodine in carbon tetrachloride simulated by Molecular Dynamics (Moldy). At the arrival of the laser flash, i.e. time zero, the background goes green. In the simulation, the atoms recombine in 40 ps after multiple collisions with the solvent molecules (shown as sticks for transparency). The time step in the movie is 0.5 ps (courtesy Rodolphe Vuilleumier, ENS, Paris). View the animation.

The experiments were carried out at the ESRF's ID09B beamline which specialises in the study of structural dynamics. The beamline is equipped with a high speed chopper that is used as a fast mechanical shutter to select the flash produced by a single bunch out of the high frequency train of flashes produced by the storage ring.

"The unique filling patterns of the storage ring, which have been developed since the outset of the institute in collaboration with the accelerator division, and specialised chopper technology, have provided us with the tools available to achieve these unique results", says M. Wulff coauthor of the publication. "We hope now to take the experiments further and observe the quantum nature of chemical bond formation, most probably in collaboration with the European X-FEL".

References

Filming the birth of molecules and accompanying solvent rearrangement, J.H. Lee, M. Wulff, S. Bratos, J. Petersen, L. Guerin, J.-C. Leicknam, M. Cammarata, Q. Kong, J. Kim, K.B. Møller and H. Ihee, *J. Am. Chem. Soc.* **135**, 3255–3261 (2013).